




PRIME
MEDICAL llc

CLINICALLY DRIVEN INNOVATION™

PRIME MEDICAL SURGICAL SOLUTIONS

Bean Bag Pressure/Protection Pad

- ▣ Provides skin protection and pressure reduction for Lateral positioning in conjunction with all types bean bag devices.
- ▣ Made of the same gel-infused viscoelastic foam found in the STP Series of Trendelenburg, reverse-Trendelenburg and Lateral positioning pads.
- ▣ Dimensions of 38" x 40" x 1" covers entire bean bag surface and allows maximum overlay coverage.
- ▣ Helps protect against nerve injury, compartment syndrome and pressure ulcers.
- ▣ Stored in rolled "sleeping bag" form for ease of storage.

BBPPP



Case Pack : One per pack

