



## PRIME MEDICAL SURGICAL SOLUTIONS

# Bean Bag Overlay

- Provides skin protection and pressure reduction for Lateral positioning in conjunction with all types of bean bag devices.
- Made of the same next-generation, high-performance polyurethane foam as found in the PrimePad Series of Anti-Skid Pads for Trendelenburg, reverse-Trendelenburg and lateral positioning.
- Dimensions of 38" x 40" x 1" covers entire bean bag surface and allows maximum overlay coverage.
- Helps protect against nerve injury, compartment syndrome and pressure ulcers.
- Stored in compressed format which lowers shipping costs and reduces storage space requirements.

Not made with natural rubber latex

PrimePad - BBO



Case Pack: 6

